

PRODUCT DESCRIPTION

The mouthpiece for the spirometer is a disposable medical device consisting of 4 layers of paper rolled into a tube shape. The outer layer of the paper is additionally coated with polyethylene in order to increase the comfort of use. The mouthpiece is used in the spirometry test, it is an element of increased hygiene of the test by protecting both the patient and the device against cross contamination.

TECHNICAL SPECIFICATIONS

<i>Material</i>	Cardboard
<i>Package</i>	Foil BBOP, cardboard box
<i>Storage</i>	Temperature 5-30°C, relative humidity <75%
<i>Product class</i>	Ila, rule 5

SIZE CHART

<i>Model</i>	<i>Device</i>	<i>Inner diameter mm.</i>	<i>Outer Diameter mm.</i>	<i>Packaging</i>
<i>PQUS-20</i>	Microlab* Bosch / Dimeq Spiro 501* SmokeCheck Micro CO Baby CO	21,6	19,6	100 pcs.
<i>PQUS-22</i>	EasyOne Schiller SpiroScout	23,9	21,9	100 pcs.
<i>PQUS-26</i>	abcMED (Pneumo) BTL Bosch / Dimeq Spiro 501	28,0	26,0	100 pcs.
<i>PQUS-28</i>	MIR (Spirolab, Spirobank, MirSpir) Vitalograph Microlab MiniWright (piklfometr) MicroPeak (pikflometr)	29,8	27,8	100 pcs.
<i>PQUS-28P</i>	abcMED (Pneumo)	30,2	27,8	100 pcs.
<i>PQUS-30</i>	MES (Lungtest)	31,0	29,0	100 pcs.
<i>PQUS-30S</i>	Aspel (Spiro 31) Jeager	31,3	29,3	100 pcs.
<i>PQUS-31</i>	abcMED (Pneumo)	31,8	29,8	100 pcs.

* with an pediatric adapter